

Test Certificate

concerning the test of influences on MRI systems caused by non-Siemens products

1. Test certificate, product designation and manufacturer

This is to certify that the following product

MR-Incubator nomag IC 3.0

model LMT-0100003-C tested with the serial number 00046,

from the manufacturer:

LMT Lammers Medical Technology GmbH

Osterweide 8

23562 Lübeck

does not affect the safety and functionality of the MRI systems of Siemens AG, Medical Solutions, specified in section 2:

Potential functional restrictions of the MRI systems are indicated in section 4.

Possible adverse effects of the MRI systems on the product indicated above are expressly not subject of the test on which this certificate is based. Therefore, this certificate does not imply non-interference of the above-mentioned product through the MRI system.

2. MAGNETOM systems involved

Cons. No.	MRI system or option (type designation)
1	MAGNETOM Trio A Tim System (3 Tesla)
2	MAGNETOM Verio (3 Tesla)

3. Applications

The nomag IC 3.0 MR-Incubator is intended for use for the following purpose in connection with magnetic resonance application; the test certificate refers to this purpose:

1. MR Diagnostics Incubator

4. Restrictions

During the use of this product the following restrictions of functionality and/or applications of the MRI systems, options and accessories apply:

1. The corresponding appropriate SpO₂- Sensor (0703007-A / LNCS YI) and SpO₂- Adapter (0703008-A / Red LNC-04) must be used only

5. Warnings


The following safety precautions must be observed when using this product in connection with the MRI systems listed in section 2:

1. The voltage transformer of the nomag IC 3.0 MR- Diagnostics- Incubator nomag IC 3.0 must be operated at a safe distance from the magnet and be equipped with a warning sign.
2. The LMT nomag TR must be placed outside the cubicle while scanning is taking place.
3. No steel cylinders may be transported with the LMT nomag TR. It and its instruction manual must bear a warning sign.
4. The material of the rubber gaiters may be visible within MR- images of MR-applications using very short spin echo times (TE << 4ms). The instruction manual must contain a corresponding warning.

6. Validity:

This test certificate shall be valid until revoked by Siemens AG, Medical Solutions.
Erlangen, 2008-06-12

Siemens AG



Weichselbaum



Schumann